


# THERALIGHT 360

## SPECIFICATIONS AND PLACEMENT

### MINIMUM ROOM DIMENSIONS




**PLEASE NOTE:** Failure to adhere to these minimum space requirements can result in compromised performance and/or damage to the TL360, thereby voiding its warranty.

### PRODUCT SPECIFICATIONS

- Length:** 87 inches (221 cm)
- Width:** 48 inches (122 cm)
- Closed Height:** 53 inches (135 cm)
- Open Height:** 68 inches (163 cm)
- Weight:** 620 lbs (281 kg)
- Weight Capacity:** 500 lbs (227 kg)
- Voltage:** 220-240 VAC 30A max 50/60 Hz
- Minimum Distance from Back Wall:** 12 inches

### DIMENSIONAL VIEW



# THERALIGHT 360

## LED, POWER AND SHIPPING SPECIFICATIONS

### LED SPECIFICATIONS



#### LED SURFACE AREA

**TOP:** 78 x 30 in. (198 x 76 cm)

**SIDE:** 78 x 14.5 in. (198 x 37 cm)

**BOTTOM:** 78 x 30 in. (198 x 76 cm)

**TOTAL LED SURFACE AREA:**

**5,811 SQ. INCHES**

#### LED DIODES

**TOTAL:** 45,000

**IRRADIANCE:** 100 mW/cm<sup>2</sup>

**LED LIFESPAN:** 100,000 hours



### POWER REQUIREMENTS

220-240 VAC 50/60HZ  
30A Single Phase  
NEMA L6-30R  
NEMA L6-30P

*PLEASE NOTE:* Failure to properly configure electrical requirements according to NEC standards can result in compromised performance and/or damage to the TL360, thereby voiding its warranty.



Plug Provided



Receptacle Required

### SHIPPING DIMENSIONS

