




DA VINCI | MEDICAL™

INTUBIX®

Quickly Insert OG Tubes


The First Time



INTUBIX® allows easy & blind guidance of an orogastric tube or any other flexible tube/probe into the gastrointestinal tract without the traditional difficulties. Once the tube or probe is in the desired position, INTUBIX is easily peeled away and discarded.

Why INTUBIX?

- ▶ Inherent but necessary design flaw in orogastric tube prevents easy traditional NG/OG placement after endotracheal intubation
- ▶ Ease and speed of placement the first time
- ▶ No need to place hand in patient's mouth
- ▶ No wasted OG tubes from misplacement
- ▶ Can be used for Dobhoff enteral tube and temperature probe placement



Features & Specs:

- ▶ Patented Smart Release™ allows Intubix to be removed while tubes/probes remain in place
- ▶ Design is rigid enough for directional control, but flexible enough to gently slide into the posterior oropharynx.
- ▶ Flanged top for quick & easy insertion of probe or tubing
- ▶ Accommodates 18 F NG/OG tubes, esophageal temperature probes and feeding tubes
- ▶ Dimensions: 13" x 5/8" DIA.
- ▶ Material: Medical grade polished silicone
- ▶ BPA & Latex Free
- ▶ One patient use. Non-sterilizable

