


# THERALIGHT FIT

## SPECIFICATIONS AND PLACEMENT

### MINIMUM ROOM DIMENSIONS



**PLEASE NOTE:** Failure to adhere to these minimum space requirements can result in compromised performance and/or damage to the TLFIT, thereby voiding its warranty.

### PRODUCT SPECIFICATIONS

**Length:** 76 inches (193 cm)

**Width:** 38 inches (96.5 cm)

**Closed Height:** 36 inches (91.5 cm)

**Open Height:** 51 inches (129.5 cm)


**Weight:** 220 lbs (100 kg)

**Weight Capacity:** 350 lbs (159 kg)

**Voltage:** 220-240 VAC 30A max 50/60 Hz

**Minimum Distance from Back Wall:** 8 inches

### DIMENSIONAL VIEW



# THERALIGHT FIT

## LED, POWER AND SHIPPING SPECIFICATIONS

### LED SPECIFICATIONS



#### LED SURFACE AREA

**TOP:** 70 x 28.5 in. (178 x 72 cm)

**BOTTOM:** 70 x 28.5 in. (178 x 72 cm)

**TOTAL LED SURFACE AREA:**

**3,990 SQ. INCHES**

#### LED DIODES

**TOTAL:** 28,000

**IRRADIANCE:** 80 mW/cm<sup>2</sup>

**LED LIFESPAN:** 50,000 hours



### POWER REQUIREMENTS

220-240 VAC 50/60HZ  
30A Single Phase  
NEMA L6-30R  
NEMA L6-30P

*PLEASE NOTE: Failure to properly configure electrical requirements according to NEC standards can result in compromised performance and/or damage to the TLFIT, thereby voiding its warranty.*



Plug Provided



Receptacle Required

### SHIPPING DIMENSIONS

