



Cryo 4S

Standing device

The CRYO 4S is one of the most advanced CRYO devices on the market today.

It combines CRYO technology and EMS muscle stimulation. Each of the 8 plates integrate cold, muscle stimulation and LED. Cryolipolysis and EMS run parallel and simultaneously in each handpiece. The Cryo 4S is not only unique in function and effect, but also in terms of quality and design. The time and intensities of both the cold temperature and EMS can be set flexibly on the ergonomic 15" touchscreen. After the cooling plates are placed on the chosen body areas and the parameters are set, the device starts CRYO / EMS body sculpting without your permanent presence. Your customer can relax and you are free for other applications in your studio.



CRYO:

- subcutaneous fat cells freeze
- fat cells are encouraged to drain using low temperatures
- white fat cells are converted into brown ones
- no suction cups required
- works with temperatures between 0 °C and 2 °C



The CRYO 4S is one of the most up-to-date cryotherapy devices currently available on the market.

It combines cryogenics with electrical muscle stimulation (EMS) and LEDs that improve the skin's appearance. Cryogenics, muscle stimulation and LEDs are integrated into each of the eight applicators, and all three treatment stimuli run at the same time and alongside each other in every handpiece.

The machine is not just unique in terms of how it works and its effects, but also in terms of its quality and design. On the ergonomic 15" touchscreen, the times and intensity of the cold and EMS can be set to whatever is appropriate for the respective customer. Once you have placed the applicators on the customer and set the parameters to suit them specifically, the machine will perform the size reduction process without you having to be present the entire time. Your customer can relax and you are free for other treatments in your studio or practice.

CRYO 4S: Innovative. Comfortable. Efficient.

As it combines three technologies, which are perfectly coordinated with each other in its eight applicators, the CRYO 4S®, our CRYO-EMS-LED, enables you to select multiple areas of the body to treat simultaneously. The CRYO 4S® can be used as the main device at your beauty salon or to extend the range of services at your beauty salon, tanning

EMS:

- muscles are encouraged to contract
- blood circulation boosted, which in turn increases the lymph flow
- fat cells are supported as they drain
- muscle compression is stimulated in a positive manner
- muscles' energy consumption is increased



salon or gym.

How does the CRYO 4S work?

- the CRYO 4S is controlled by an intelligent piece of software that we developed ourselves in Germany
- you select the cryogenic-, EMS- and LED-parameters you want to use on the 15" touchscreen
- it can be set to suit each customer specifically using various EMS and cryogenic programmes
- treatments can involve one to eight handpieces
- different areas of the body can be treated simultaneously
- it provides safer treatments without an operator, as a PIN code prevents unauthorised changes

Our service at purchase?

- own production and workshop in Germany
- comprehensive advice
- trainings including certificate
- professional maintenance and repair
- different financing methods with beneficial conditions
- 2 years of warranty



LED:

- LED light in the 400 nm range
- improvements to the skin's appearance



MADE in GERMANY

We have developed the CRYO 4S ourselves and designed it to meet the needs of tanning salons, practices and beauty salons perfectly. In addition, all our devices are completely manufactured in Germany, guaranteeing you quality you can rely on.

Place your trust in our many years of experience and our expertise in body sculpting using cryogenics and EMS.

For more information about our Cryo 4S device or a non-binding offer, please feel free to contact us via telephone or email.

We are looking forward to hearing from you soon.

